

## MST mini



Magnetic stirrer without heating

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	0.8
Motor rating input [W]	4	Motor rating output [W]	0.5
Speed range [1/min]	0 - 1500	Speed control [-]	stepless
Set-up plate material [-]	polyester	Set-up plate dimensions [mm]	115 x 115
Dimensions (W x H x D) [mm]	115 x 40 x 130	Weight [kg]	0.23
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 42	Voltage [V]	100 - 240
Frequency [Hz]	50/60	Power input [W]	4