

## yellow MAG MS 7



New magnetic stirrer without heating

- Ceramic top plate
- Excellent resistance to chemicals
- Microprocessor technology: guarantees maximum speed control accuracy

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H2O) [l]	10
Stirring quantity max. (H2O) [l]	10	Motor rating input [W]	15
Motor rating output [W]	1.5	Speed display [-]	scale
Speed range [1/min]	100 - 1500	Stirring bar length max. [mm]	80
Speed control [-]	scale 0 - 6	Set-up plate material [-]	ceramic
Set-up plate dimensions [mm]	180 x 180	Dimensions (W x H x D) [mm]	215 x 105 x 330
Weight [kg]	5	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 21
Voltage [V]	230 / 120 / 100	Frequency [Hz]	50/60
Power input [W]	30		