

MSH basic



Magnetic stirrer with stainless steel heating plate.

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H2O) [l]	5
Stirring quantity max. (H2O) [l]	5	Motor rating input [W]	15
Motor rating output [W]	2	Speed range [1/min]	0 - 2000
Heat output [W]	400	Heating temperature range [°C]	room temp. - 350
Heat control [-]	stepless	Speed control [-]	stepless
Set-up plate material [-]	stainless steel 1.4301	Set-up plate dimensions [mm]	Ø 125
Dimensions (W x H x D) [mm]	160 x 105 x 210	Weight [kg]	2.3
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 21	Voltage [V]	230 / 115 / 100
Frequency [Hz]	50/60	Power input [W]	420